

Full-electric/Semi-electric bed

adjustable height



- Durable split spring construction
- Top grade helical springs and wire fabric
- Scratch resistant polymer epoxy coated frame
- Synchronized motor system with automatic shut-off memory
- Durable 3" self lubricating ball bearing casters, two with wheel locks
- Bed utilizes 35" x 80" x 7" mattress
- Optional bed rail and mattress

	min. bed height (in.)	max. bed height (in.)	overall length (in.)	overall width (in.)	weight capacity (lb.)
860A	16	24.5	88	35	400
860B	16	24.5	88	35	400

Model No.	Description
860A	Full-electric bed, adjustable height
860B	Semi-electric bed, adjustable height