



## Patient lift

Electric

- LINAK actuator
- Electronic controller box
- Aluminum hydraulic pump with quick release valve
- 2" square powder coating steel tubing
- 4" swivel casters with two brakes on back casters
- Foldable mast and base
- Adjustable U-shaped base
- Push bar

	height adjustment (in.)	base adjustment (in.)	overall width (in.)	total weight (lb.)	weight capacity (lb.)
895	27 to 72	24 to 38	24	94.6	330

Model No.	Description
895	Electric patient lift