



Freedom

companion wheelchair series

- Reinforced chrome-plated or powder-coated steel frame with four 8" caster
- High-grade breathable upholstery, black, padded
- Fold-down hinged back with push-button release pin for easy transport
- Padded armrests and calf supports
- Chrome-plated, metallic burgundy or blue finish

	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall width open (in.)	weight w/o rigging (lb.)	weight w/o rigging & wheels (lb.)	weight capacity (lb.)
100	17	16	18.5	18.5	39	21.5	24	n/a	220
100L	19	16	18.5	18.5	39	23.5	24	n/a	220
100L/SP	19	16	18.5	18.5	39	23.5	25	n/a	220

Model No.	Description
100	17", fixed arms, swingaway footrests
100E	17", fixed arms, elevating legrests
100L	19", fixed arms, swingaway footrests
100LE	19", fixed arms, elevating legrests
100L/SP	19", fixed arms, swingaway footrests, with side panel
100LE/SP	19", fixed arms, elevating legrests, with side panel