



Explorer

extra wide, heavy duty wheelchair series

- Heavy-duty, durable chrome-plated steel frame, double crossbar
- Reinforced rear wheel axle & bushing
- Superior fork with premium bearing, cone design similar to mountain bike head set
- 24" mag wheels with high performance flat-free tires and handrims
- 8" x 2" caster wheels with PU tires or 8" x 1" caster wheels with rubber tires
- Flame retardant double embossed upholstery, dark blue
- Striped side panels and anti-slip rubber footplates
- Padded armrests and calf supports

	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall width open (in.)	weight w/o rigging (lb.)	weight w/o rigging & wheels (lb.)	weight capacity (lb.)
130XL	22	17	20	16.5	36	30.5	50.5	n/a	400
130XXL	24	17	20	16.5	36	32.5	51.5	n/a	400

Model No.	Description
130XL	22", detachable arms, swingaway footrests
130XLE	22", detachable arms, elevating legrests
130XXL	24", detachable arms, swingaway footrests
130XXLE	24", detachable arms, elevating legrests