



## 910 - Deluxe Pump & Pad System

by American Bantex

### **Deluxe Pump (Adjustable pressure control)**

The pump offers high capacity output and features an adjustable pressure control which gives variable pressure to accommodate the patient's weight. It includes an illuminated on/off switch and a hook for attachment to the bed.

- Power: 115V/60Hz or 230V/50Hz
- Pressure Range: 50 mmhg - 105 mmhg
- Air Output: 4 - 5 liters (115V) / 5 - 6 liters (230V)
- Cycle Time: 5 minutes (60 Hz) / 6 minutes (50 Hz)
- Power Cord: 11 feet     ■ Size: 10 x 5-1/8 x 3-7/8"     ■ Weight: 3.5 lbs

### **Deluxe Bubble Pad**

The pad is made of heavy duty, bubble textured PVC that can tuck underneath a mattress at both ends. The heavy gauge premium beige vinyl resists tearing and punctures.

- Size: 118" x 34"     ■ Height of Cell: 2-1/2"
- Number of Cells: 130     ■ Thickness: 11 gauge