



920 - Therapy Pump & Low Air Loss Mattress System

by American Bantex

Therapy Pump

- Power: 115V/60Hz or 230V/50Hz
- Pressure Range: 40 mmhg - 90 mmhg
- Air Output: 4 - 5 liters (115V) / 5 - 6 liters (230V)
- Cycle Time: 8 minutes (60 Hz) / 10 minutes (50 Hz)
- Power Cord: 11 feet ■ Size: 10 x 5-1/8 x 3-7/8" ■ Low pressure light

Low Air Loss Mattress

- Size: 80" x 36" x 5" ■ Air Cell Depth: 5"
- 17 Nylon and TPU air cells ■ Nylon PU top cover
- CPR valve for quick emergency air release
- 8 air cells with laser holes for air loss
- Removable non-shearing quilted waterproof top cover encapsulates the entire mattress
- Replacement air cells available