

Bravo Series

shroud model, power wheelchair



- Front wheel drive shroud model power chair
- 18" reclining, high back seat
- Powerful 450 W motors
- The arms flip back and adjust in and out, for maximum user comfort
- Adjustable seat height
- The footplate flips up and is height adjustable
- Safety seat belt and charger included; battery excluded
- Candy apple red

| | | | |
|--------------------------|-------------------------------|---|--------------------------------------|
| DIMENSION | Length: 40.2" Width: 24.4" | SPEED RANGE TURNING RADIUS | 4.5 mph 25 miles 24.2" |
| SEAT HEIGHT | 19.0" | INCLINE CAPABILITY | 12° |
| | | GROUND CLEARANCE | 2.8" |
| DRIVE WHEELS | 10" X 3.3" Solid | BRAKE | Full-Auto Magnetic Disc Safety Brake |
| REAR WHEELS | 8" X 2 Solid | WEIGHT W/BATTERY | 203.5 lbs |
| CARRYING CAPACITY | 300 lbs | WEIGHT W/O BATTERY | 138.4lbs |

| Model No. | Description |
|-----------|---|
| HP4F/BURG | 18" mid wheel drive powerchair, flip-up footplate |

