

# Bravo Series

shroud model, power wheelchair



- Compact mid-wheel-drive power base
- Dynamic shark controller with extendable/retractable remote joystick
- 18" mid-back swivel seat
- The arms flip back and adjust in and out, for maximum user comfort
- Easily accessible free-wheeling lever
- Built in charger
- The footplate flips up and is height adjustable
- Safety seat belt and charger included; battery excluded
- Candy apple red, silver & blue

<b>DIMENSION</b>	Length: 37.5" Width: 23.3" Height: 40.3"	<b>SPEED RANGE</b> <b>TURNING RADIUS</b>	4.5 mph 25 miles 30.7"
<b>SEAT HEIGHT</b>	20 - 1/2" to 22 - 1/2"	<b>INCLINE CAPABILITY</b>	8°
<b>DRIVE WHEELS</b>	10" X 2 Solid	<b>GROUND CLEARANCE</b>	2.4"
<b>REAR WHEELS</b>	8" X 2 Solid	<b>DRIVE TRAIN</b>	2-Motor Mid-wheel drive, 420W
<b>WEIGHT W/O BATTERY</b>	120 lbs	<b>CARRYING CAPACITY</b>	300 lbs

Model No.	Description
HS1000	19" rear wheel drive powerchair, flip-up footplate



Medicare HCPCS: K0011

