



Nautilus

four wheel rollator

- Heavy-duty durable steel frame
- Adjustable height handles
- Anatomic handgrips
- 8" x 2" EVA tires
- Large seat (15.5" x 9.0")
- Removable wire basket, plastic tray included
- Locking loop
- Metallic blue finish

	seat height (in.)	seat size (in.)	handle height (in.)	width between handles (in.)	overall width (in.)	width open (in.)	width closed (in.)	net weight (lb.)	weight capacity (lb.)
515	24.0	15.5 x 9.0	30.0 - 38.0	16.5	23.0	-	-	22.5	250

Model No.	Description
515	locking loop brakes