



# Trident

three wheel walker

- Durable steel frame
- Adjustable height handles
- Anatomic handgrips
- 8" x 2" EVA tires
- Locking loop brake
- Plastic bag included
- Basket optional
- Chrome-plated, metallic burgundy or blue finish

	seat height (in.)	seat size (in.)	handle height (in.)	width between handles (in.)	overall width (in.)	width open (in.)	width closed (in.)	net weight (lb.)	weight capacity (lb.)
510	-	-	31.0 - 39.0	17.0	-	27.0	10.5	11.9	225

Model No.	Description
510	locking loop brakes