

Snazzy Series

3 wheel scooter, economic



- The programmable control panel has easy-to-read battery gauge and self-diagnosing flash unit
- Scooter can be dismantled quickly i.e. no tools required
- Swivel seat with adjustable height can be folded easily
- Thumb operates forward & reverse control, and speed adjustable knob
- 9" front and rear driving wheels with flat-free tires offer excellent maneuverability even in tight spaces
- The portable on-board charger conveniently stores in the trunk compartment and goes everywhere with the chair
- Reflectors and anti-tipper improve safety
- Optional mirror, oxygen tank holder and cane holder
- Candy apple red, emerald green or silver

| | | | |
|--------------------------|--|---------------------------|---------------------------|
| DIMENSION | Length: 43.0" Width: 21.0" Height: 37.0" | SPEED RANGE | 5 mph 15 miles 3.4" |
| CARRYING CAPACITY | 250 lbs | GRADE CLIMBABLE | 12° |
| WHEELS | 9" diameter X 3 | CURB CLIMBING | 3.2" |
| BRAKING SYSTEM | Double automatic safety braking systems | PEAK MOTOR POWER | 560W 3/4HP |
| BATTERY CHARGER | 24V, 1.8A | WEIGHT W/O BATTERY | 100.8 lbs |

| Model No. | Description |
|-----------|---|
| HS320 | 18", three-wheel power mobility scooter, economic |



Medicare HCPCS: E1230

