

Snazzy Series

3 wheel scooter, economic



- Maximum Speed: 4 mph
- Capacity: 250 lbs
- The 18.7" wide track provides more stability when negotiating corners.
- Stronger, reinforce chassis to reduce flex and rolling.
- World class electronic controls ensure safety and peace of mind for the riders.
- One-Piece battery pack module for ease of disassembly.
- Padded foldable seat and expandable legroom to provide a custom fit.
- Standard 12Ah/10mile range batteries
- Yellow and Red color

DIMENSION	Length: 37.6" - 40.0" Width: 18.7" Height: 33.8"	SPEED RANGE GROUND CLEARANCE	4 mph 10 miles 1.6"
CARRYING CAPACITY	250 lbs	GRADE CLIMBABLE	8°
WHEELS	8" diameter X 3	CURB CLIMBING	1.9"
BRAKING SYSTEM	Electro-Mechanical	MOTOR SIZE	250W, 4000 r.p.m.
BATTERY CHARGER	2A Built In	WEIGHT W/O BATTERY	63 lbs
SUSPENSION	N/A	TURNING RADIUS	35.4"

Model No.	Description
HS120	14", three-wheel power mobility scooter, economic



Medicare HCPCS: E1230

