

# Snazzy Series

3 wheel scooter, economic



- Maximum Speed: 4 mph
- Capacity: 250 lbs
- The 18.7" wide track provides more stability when negotiating corners.
- Stronger, reinforce chassis to reduce flex and rolling.
- World class electronic controls ensure safety and peace of mind for the riders.
- One-Piece battery pack module for ease of disassembly.
- Padded foldable seat and expandable legroom to provide a custom fit.
- Standard 12Ah/10mile range batteries
- Yellow, Blue and Red color

<b>DIMENSION</b>	Length: 39.1" - 41.5" Width: 18.7" Height: 33"	<b>SPEED RANGE</b> <b>GROUND CLEARANCE</b>	4 mph 10 miles 1.6"
<b>CARRYING CAPACITY</b>	250 lbs	<b>GRADE CLIMBABLE</b>	8°
<b>WHEELS</b>	8" diameter X 3	<b>CURB CLIMBING</b>	2"
<b>BRAKING SYSTEM</b>	Electro-Mechanical	<b>MOTOR SIZE</b>	250W, 4000 r.p.m.
<b>BATTERY CHARGER</b>	2A Built In	<b>WEIGHT W/O BATTERY</b>	66.7 lbs
<b>SUSPENSION</b>	N/A	<b>TURNING RADIUS</b>	35.8"

Model No.	Description
HS235	14", three-wheel power mobility scooter, economic



Medicare HCPCS: E1230

