

Snazzy Series

4-wheel scooter



- The programmable control panel has easy-to-read battery gauge and self-diagnosing flash unit
- Scooter can be dismantled quickly i.e. no tools required
- Swivel seat with adjustable height can be folded easily
- Thumb operates forward & reverse control, and speed adjustable knob
- 10" front and rear driving wheels with flat-free tires offer excellent maneuverability even in tight spaces
- The portable on-board charger conveniently stores in the trunk compartment and goes everywhere with the chair
- The front and rear bumper prevents the scooter from bumping damage
- Headlight, rear reflectors, and anti-tipper improve safety
- Optional mirror, oxygen tank holder and cane holder
- Candy apple red, emerald green or silver

DIMENSION	Length: 50.0" Width: 24.0" Height: 39.0"	SPEED RANGE GROUND CLEARANCE	5 mph 25 miles 3.0"
CARRYING CAPACITY	300 lbs	GRADE CLIMBABLE	12° ± 1°
WHEELS	10" diameter X 4	TURNING RADIUS	62.6"
BRAKING SYSTEM	Double automatic safety braking systems	MOTOR	DC 24V, 370 watt, 3400 rpm
BATTERY CHARGER	24V, 4AH	WEIGHT W/O BATTERY	140.4 lbs

Model No.	Description
HS580	18", four-wheel power mobility scooter



Medicare HCPCS: E1230

